

## RSAR001 user's manual

### Culture Media

- RPMI1640 medium with 15% heat-inactivated FBS

### Reagents

- Sterilized D-PBS
- PBS-EDTA
- Trypsin-EDTA\*
- Trypsin neutralizing solution\*
- HEPES buffered saline solution\*

\*Recommend to use "ReagentPack™ Subculture Reagents" (Lonza, CC-5034)

- Cryopreservation medium

### Handling procedure for frozen cells

1. Add 15 ml of medium into 50-ml sterile tube in a biological safety cabinet.
2. Thaw the cells by gentle agitation in a 37 °C water bath for 1 minute.
3. Take the vial out of the water bath, wipe with 70% ethanol, then into the biological safety cabinet.
4. Transfer the vial contents into a sterile tube containing 15 ml of medium.
5. Mix the cells and medium gently.
6. Transfer the cell suspension to a T75 flask.
7. Incubate the culture at 37 °C, 5% CO<sub>2</sub>.
8. Next day, remove all culture medium from the T75 flask, and add 15 ml of fresh medium.

### Subculturing procedure

1. Remove culture medium from the T75 flask.
2. Rinse twice the cell layer with 10 ml of PBS-EDTA.
3. Rinse the cell layer with 6 ml of HEPES buffered saline solution.
4. Add 5 ml of Trypsin-EDTA to flask and incubate at room temperature for 5 minutes.
5. Add 10 ml of Trypsin neutralizing solution and gently pipette the cell suspension.
6. Transfer the culture suspension to a sterile centrifuge tube and centrifuge at 200 x g for 3 minutes at room temperature.
7. Discard the supernatant.
8. Resuspend the cell pellet in 10 ml of D-PBS.
9. Centrifuge at 200 x g for 3 minutes at room temperature.

10. During centrifugation, add 14 ml of medium into a T75 flask.
11. Discard the supernatant.
12. Resuspend the cell pellet in 3 ml of fresh medium and add 1 ml of the cell suspension into the T75 flask containing 14 ml of medium (1/3 split).
13. Incubate culture at 37 °C.

#### Cryopreservation procedure

1. Follow the step 1 to 9 of “Subculturing procedure”.
2. Discard the supernatant.
3. Resuspend the cell pellet in 3 ml of cryopreservation medium.
4. Add 1 ml each of cell suspension into three cryopreservation vials.
5. Place the vial in a cryopreservation container and freeze in a -80 °C freezer.
6. Store the cryopreserved vial in vapor phase of liquid nitrogen.