CHO DC-SIGN-N80A cell line

Catalogue number: 160717

Sub-type: Images:

Contributor

Inventor: Alessandra Cambi

Institute: Radboud University Medical Center

Images:

Tool details

*FOR RESEARCH USE ONLY

Name: CHO DC-SIGN-N80A cell line

Alternate name: DC-SIGN N8A

Class: Conjugate:

Cancer Tools.org Description: Dendritic Cell-Specific ICAM-3-Grabbing Non-Integrin (DC-SIGN; CD209) is a type II plasma membrane PRR abundantly expressed in antigen-presenting cells such as dendritic cells (DCs) and activated macrophages. As member of the C-type Lectin Receptor (CLR) family, DC-SIGN binds a plethora of pathogens, ranging from viruses like HIV-1, ebola virus, and hepatitis C virus, to larger pathogens like Mycobacterium tuberculosis and Candida albicans. This N80A mutant lacks the Nglycosylation ligand binding domain.

Purpose:

Parental cell: CHO Organism: Hamster

Tissue: Ovary

Model: Gender: Isotype: Reactivity: Selectivity: Host:

Immunogen:

Immunogen UNIPROT ID:

Sequence:

Growth properties:

Production details: CHO cell lines stably expressing DC-SIGN-N80 were established by Lipofectamin

2000 transfection, and were cultured in HamÄ?Ë???Â???Â?s F-12 medium supplemented with 10% heat-inactivated FBS, 1% Antibiotic Antimycotic Solution, and 0.5Ä?Ë???Â?mg/ml of the aminoglycoside antibiotic G418.

Formulation:

Recommended controls: Bacterial resistance: Selectable markers: Additional notes:

Target details

Target: Dendritic Cell-Specific Intercellular adhesion molecule-3-Grabbing Non-integrin lacking the ligand binding domain N80.

Cancer Tools.org

Target alternate names:

Target background:

Molecular weight:

Ic50:

Applications

Application:

Application notes:

Handling

Format: Frozen
Concentration:
Passage number:

Growth medium: HamÄ?Ë???Â???Â?s F-12 medium supplemented with 10% heat-inactivated FBS, 1% Antibiotic Antimycotic Solution , and 0.5Ä?Ë???Â?mg/ml of the aminoglycoside antibiotic

G418.

Temperature: Atmosphere:

Volume:

Storage medium: Storage buffer:

Storage conditions:

Shipping conditions: Dry ice

Related tools

Related tools: CHO DC-SIGN WT cell line ; CHO DC-SIGN-Ä?Â???Â?Rep cell line ; CHO DC-SIGN-Ä?Â???Â?CRD cell line

References

References:

