HeLa GRASP55 knockout cell line

Catalogue number: 157706 Sub-type: Continuous

Images:

Contributor

Inventor: Yanzhuang Wang Institute: University of Michigan

Images:

Tool details

*FOR RESEARCH USE ONLY

Name: HeLa GRASP55 knockout cell line

Alternate name:

Class:

Conjugate:

Cancer Tools.org **Description:** Clustered regularly interspaced short palindromic repeats (CRISPR)/Cas9 technology was used to knock out Golgi reassembly stacking protein of 65 kDa (GRASP65) and/or Golgi reassembly stacking protein of 55 kDa (GRASP55) in HeLa and HEK293 cells. By knocking out GRASPs it was demonstrated that these proteins are essential for Golgi stacking structure formation and accurate posttranslational modifications function in the Golgi. The GRASP knock out cell lines can be used to study the role of GRASPs in Golgi-related biological processes and diseases in which the Golgi is defective (e.g. Alzheimers disease, congenital disorders of glycosylation, foot-and-mouth disease, reoxygenation injury and cancer). The HeLa WT cell line can be provided, as a positive control, if requested with order.CRISPR edited HeLa cells.

Purpose:

Parental cell: HeLa cell line

Organism: Human Tissue: Cervix Model: Knock-Out

Gender: Isotype: Reactivity: Selectivity: Host:

Immunogen:

Immunogen UNIPROT ID:

Sequence:

Growth properties:

Production details: Duplexed sgRNA oligos were digested and ligated into pSpCas9(BB)-2A-GFP(PX458) and pSpCas9(BB)-2A-Puro(PX459)v plasmids to generate GRASP55/65 GFP or Puro plasmids, respectively. CRISPR knock out cells were generated by transfection with GRASP55/65 GFP or Puro plasmids followed by enrichment of GFP-expressing cells by FACS sorting or by selection with 1 Ä?Â?????g/ml puromycin, respectively. Individual clones were generated by plating cells at low density and isolating individual colonies. GRA...

Formulation:

Recommended controls: HeLa parental line

Bacterial resistance: Selectable markers:

Additional notes: The HeLa WT cell line can be provided, as a positive control, if requested with order. CRISPR edited HeLa cells. Cancer Research Technology Limited (trading research tools as Ximbio) has been granted a non-exclusive license to the CRISPR-Cas9 technology by ERS Genomics Ltd under the patent rights listed here. This license from ERS Genomics Ltd allows Ximbio to develop an...

Cancer Tools.org

Target details

Target: GRASP55

Target alternate names:

Target background:

Molecular weight:

Ic50:

Applications

Application:

Application notes: Cancer Research Technology Limited (trading research tools as CancerTools.org) has been granted a non-exclusive license to the CRISPR-Cas9 technology by ERS Genomics Ltd under the patent rights listed here: https://www.cancertools.org/tool-

faqs#hs_cos_wrapper_widget_1649861453796 This license from ERS Genomics Ltd allows CancerTools.org to develop and commercialise CRISPR-Cas9 modified cell lines for research use only. CancerTools.org can provide these modified CRISPR-Cas9 cell lines to comp...

Handling

Format: Frozen
Concentration:
Passage number:
Growth medium:

Temperature: Atmosphere:

Volume:

Storage medium: Storage buffer:

Storage conditions:

Shipping conditions: Dry ice

Related tools

Related tools:

References

References: Bekier ME et al. 2017. Mol Biol Cell. 28(21):2833-2842. PMID: 28814501.